

# IN CONTINENCE AND UROLOGY PRODUCTS

# FOLEY CATHETERS



## Description

Size : 6,8,10Fr/Ch  
 Balloon Capacity : 6Fr/Ch 3cc,8,10Fr/Ch 3-5cc  
 Length : 250mm  
 Made from natural latex. Silicone Coated.  
 Available with the option of plastic valve or rubber valve.  
 Natural or Hydrophilic.  
 Packed sterile. For Single use only.



## Description

Size : 12,14,16,18,20,22,24,26,28,30Fr/Ch  
 Balloon Capacity : 12-24Fr/Ch 5-10cc,5-15cc  
 : 12-30Fr/Ch 30cc  
 : 16-30Fr/Ch 30-50cc  
 Length: 400mm(Male) or 270mm(Female)  
 Made from natural latex. Silicone Coated.  
 Natural or Hydrophilic and Male or Female.  
 Packed sterile. For Single use only.



## Description

Size : 16,18,20,22,24,26Fr/Ch  
 Balloon Capacity : 30cc  
 Length : 400mm  
 Made from natural latex. Silicone Coated.  
 Natural or Hydrophilic.  
 Packed sterile. For Single use only.

# SILICONE FOLEY CATHETERS



## Description


Size : 8,10,12,14,16,18,20,22,24Fr/Ch  
 Balloon Capacity : 8-10Fr/Ch 3cc,3-5cc,5cc  
 : 12-24Fr/Ch 5-10cc,5-15cc,5cc,10cc,15cc  
 : 16-24Fr/Ch 30cc  
 Length : 400mm  
 Made from 100% Silicone.  
 Packed sterile. For Single use only




## Description

Size : 16,18,20,22,24Fr/Ch  
 Balloon Capacity : 5-15cc,30cc  
 Length : 400mm  
 Made from 100% Silicone.  
 Packed sterile. For Single use only.

# SUCTION CATHETER



# NELATON CATHETER



## Description

Size : 6,8,10,12,14,16,18,20,22,24Fr/Ch  
Length : 400mm(Male) or 180mm(Female)  
Made from PVC or SOFT PVC.  
Natural or Hydrophilic and Male or Female.  
Packed sterile. For Single use only.

# ECONOMY URINE BAG

## Description



Ref : GCU201501  
Front Film : Transparence  
Back Film : White  
Outlet : Pull-push valve  
Capacity : 2000ml  
A. 1500ml and 1000ml are available.  
B. With anti-reflux valve. Tubing length will be 90cm or 110cm.  
C. Available in two plastic hangers or belts. Different distance is available between double holes.  
D. Latex-free. Packed in poly bag or blister with E.O. sterile. Or bulk. Greatcare brand or OEM available.



Ref : GCU201505 With reinforced plastic hang tape  
Front Film : Transparence  
Back Film : White  
Outlet : Pull-push Valve  
Capacity : 2000ml



Ref : GCU201701  
Front Film : Transparence  
Back Film : White  
Outlet : Without  
Anti-reflux valve: Without  
Capacity : 2000ml

Ref : GCU201703  
Front Film : Transparence  
Back Film : White  
Outlet : Without  
Anti-reflux valve: With  
Capacity : 2000ml

A. 1500ml and 1000ml are available.  
B. With anti-reflux valve. Tubing length will be 90cm or 110cm.  
C. Available in two plastic hangers or belts. Different distance is available between double holes.  
D. Latex-free. Packed in poly bag or blister with E.O. sterile.  
Or ulk. Greatcare brand or OEM available.

# DELUXE URINE BAG

## Description

Urine Bag  
2000ml,with white sample port & T valve.  
Extra thickness film(0.15mm per ply),double transparent films.  
Extra outer diameter tubing(9mm).  
Four holes for different distance hang sizes.

## Description

Urinary Drain Bag  
2000ml,with T valve.  
Front film is transparent, back film is white.  
With white sample port.With anti-reflux valve.  
Latex-free.Packed sterile.

## Description

Urinary Drain Bag  
2000ml,with twist valve.  
Front film is transparent, back film is white.  
With sample port. With anti-reflux valve. With plastic tie.  
Latex-free.Packed sterile.

## Description

Urinary Drain Bag  
2000ml,with T valve.  
Front film is transparent, back film is white.  
With sample port.With anti-reflux valve.With plastic tie.  
Latex-free.Packed sterile.

## Description

Urinary Drain Bag  
2000ml,with twist valve.  
Front film is transparent, back film is white.  
With sample port.With anti-reflux valve. With reinforced single hanger.  
Latex-free.Packed sterile.





# DELUXE URINE BAG



## Description

Urinary Drain Bag  
2000ml, with twist valve.  
Front and back film are transparent.  
With anti-reflux valve. With plastic tie.  
Latex-free. Packed sterile.



## Description

Urinary Drain Bag  
2000ml, with T valve. With white sample port.  
Front film is transparent, and back film is white.  
With anti-reflux valve. With reinforced double hanger. With tubing clamp. Elastic tubing.  
Latex-free. Packed sterile



## Description

Urinary Drain Bag  
Ref.: GCU210701 2600ml  
Ref.: GCU210702 2500ml  
Ref.: GCU210703 2000ml  
With anti-reflux tower (three part), With needless sample port.  
With bed sheet clamp. Tubing clip is optional.  
With reinforced double hanger & rope hanger.  
With T valve. Latex-free. Packed sterile.  
Blister pack or poly bag packing.




## Description

Urinary Drain Bag  
Ref.: GCU210801 2600ml  
Ref.: GCU210802 2500ml  
Ref.: GCU210803 2000ml  
With anti-reflux tower (one part), With needless sample port.  
With bed sheet clamp. Tubing clip is optional.  
With reinforced double hanger & rope hanger.  
With T valve. Latex-free. Packed sterile.  
Blister pack or poly bag packing.



## Description

Urinary Drain Bag  
2000ml, With anti-reflux tower, With needless sample port.  
With reinforced single hanger. With finger guard clamp.  
Latex-free. Packed sterile.  
Blister pack or poly bag packing.



# DELUXE URINE BAG

## Description



Urinary Drain Bag  
2000ml, With anti-reflux tower, With white sample port.  
With reinforced single hanger. With T valve outlet.  
Latex-free.Packed sterile.  
Blister pack or poly bag packing.

## Description



Urinary Drain Bag  
Ref: GCU211701 2600ml  
Ref: GCU211702 2500ml  
Ref: GCU211703 2000ml  
With anti-reflux device(two parts), with needless sample port.  
With bed sheet clamp. Tubing clip is optional.  
With reinforced double hanger & rope hanger.  
With T valve.Latex-free.Packed sterile.  
Blister pack or poly bag packing.

# URINE METER



Ref.: GCU220101 2600ml drain bag and 200ml urine meter.  
Ref.: GCU220102 2500ml drain bag and 200ml urine meter.  
Ref.: GCU220103 2600ml drain bag and 400ml urine meter.  
A. Easy-to -read molded plastic meters are calibrated up to 200ml.  
B. With anti-reflux valve.With needless sample port.  
C. With reinforced double hanger and rope hanger.  
D. Latex-free.Packed sterile.



2600ml drain bag and 500ml urine meter.



2600ml drain bag and 400ml urine meter.  
With big reinforced double hanger & rope hange.



## DELUXE URINE BAG




2600ml drain bag and 500ml urine meter.  
With drip chamber and white sample port.



2000ml drain bag and 200ml urine meter.

## LEG BAG



500ml, ellipse shape, with  
green twist valve &  
green tip cover.  
With double high elastic belts.



600ml, square shape, with  
green pull-push valve &  
green tip cover. With double high  
elastic belts.



500ml, one ply flocking fabric,  
another ply is PVC film.  
With T valve & short tubing.




500ml, one ply flocking fabric,  
another ply is PVC film.  
With T valve & long tubing.




500ml, one ply flocking fabric,  
another ply is PVC film.  
With double reversal valve &  
long tubing.

## URINE COLLECTOR



Pediatric Urine Collector  
100ml, frankseal,  
yellow paper.



Pediatric Urine Collector  
100ml, frontseal, yellow  
paper.



Pediatric Urine Collector  
160ml, frankseal, white  
paper.



Enema Set 1500ml



## WOUND DRAINAGE RESERVOIR



Wound Drainage Reservoir  
With anti-reflux valve, with different sizes needle,  
such as 8fr, 10fr, 12fr, 14fr, 16fr, 18fr, etc.





# URETHRAL CATHETER



Size: ,8,10,12,14,16,18,20,22,24,26,28Fr/Ch  
Length: 400mm  
Packed sterile.For Single use only.



Size: 6,8,10,12,14,16,18,20,22,24,26,28Fr/Ch  
Length: 400mm  
Packed sterile.For Single use only.

## MALECOT CATHETER



Size: 12,14,16,18,20,22,24,26,28,30,  
32,34,36,38,40Fr/Ch  
Length: 400mm  
Packed sterile.For Single use only.

## MALE EXTERNAL CATHETER



Packed sterile, For single use only



100% silicone,self-adhering,  
with white tray.  
transparent tray available



100% Silicone,Self-adhering  
Type: Standard,Pop on, Wide Band  
Size: 24mm,28mm,31mm,35mm,40mm

## URINE BAG HANGER



## URINE BAG VALVE



## URINE BAG CLAMP



## URINE BAG CONNECTOR



## DIAPER



DIAPER X-LARGE WITH INDICATOR  
DIAPER LARGE WITH INDICATOR, 80/CTN  
DIAPER SMALL WITH INDICATOR, 80/CTN  
DIAPER MEDIUM WITH INDICATOR, 80/CTN  
DIAPER (PAMPERS) 16+  
DIAPER (PAMPERS) 9-20KG  
DIAPER (PAMPERS) 4-9 KG  
DIAPER (PAMPERS) 2-5 KG  
DIAPER NEW BORN 3-6 KG 64 X 3/CTN  
DIAPER(PAMPERS) 6X18/CTN

## UNDERPAD



UNDER PAD 60 X 45 CM 38G  
UNDER PAD 60 X 60 CM 60G-240/CTN  
UNDER PAD 60 X 90 CM70G-150/CTN  
UNDER PAD 60 X 90 CM80G-150/CTN





## OXYGEN AND AIROSOL THERAPY OXYGEN MASK



### Characteristics

OXYGEN MASK TRANSPARENT-ADT  
OXYGEN MASK TRANSPARENT-CHILD

## NEBULIZER MASK



### Characteristics

NEBULIZER MASK INFANT  
NEBULIZER MASK ADULT  
NEBULIZER MASK CHILD  
NEBULIZER MASK ADULT TRIDENT  
NEBULIZER MASK CHILD TRIDENT

## VENTIMASK



## VENTURI MASK



### Characteristics

VENTURI MASK ADULT TRANSPARENT  
VENTURI MASK CHILD GREEN

## NON REBREATHER MASK



### Characteristics

NON REBREATHER MASK TR.ADULT  
NON REBREATHER MASK TRA.CHLD

# ANESTHESIA & RESPIRATORY

## ANESTHESIA MASK



### Characteristics

ANESTHESIA MASK ADULT  
ANESTHESIA MASK ASSORTED  
ANESTHESIA MASK CHILD  
ANESTHESIA MASK-SIZE 1  
ANESTHESIA MASK-SIZE 2  
ANESTHESIA MASK-SIZE 3  
ANESTHESIA MASK-SIZE 4  
ANESTHESIA MASK-SIZE 5  
ANESTHESIA MASK-SIZE 6

## OXYGEN T-PIECE



## TRACHOSTOMY MASK




## VENTIMASK



## HUMIDITY ADAPTOR



## T-PIECES, TUBING WATER TRAPS & CONNECTORS



## NASAL CANNULA



### Characteristics

NASAL CANULA ADLT TRANSPARENT  
NASAL CANULA INFANT TRANSPARENT  
NASAL CANULA CHLD TRANSPARENT

# ANESTHESIA & RESPIRATORY


## AMBU BAG & MASK



### Characteristics

AMBU MASK 1-ADLT LGE 1PIECE AIR CUS  
 AMBU MASK 3-CHILD 1PIECE AIR CUS  
 AMBU MASK 4-INFANT 1 PIECE AIR CUS  
 AMBU MASK-5 ADLT 1 PIECE AIR CUS  
 AMBU BAG PVC DISPOSABLE ADULT  
 AMBU BAG PVC DISPOSABLE CHILD  
 AMBU BAG PVC DISPOSABLE INFANT  
 AMBU BAG SILICONE RESUSCITATOR-ADULT  
 AMBU BAG SILICONE RESUSCITATOR-CHILD  
 AMBU BAG SILICONE RESUSCITATOR-INFANT

## INTUBATION STYLET



## INTUBATION STYLET



## NASOPHARYNGEAL AIRWAY



GUEDEL AIRWAYS - SIZE #0-BLACK  
 GUEDEL AIRWAYS - SIZE #00  
 GUEDEL AIRWAYS - SIZE #000-PINK  
 GUEDEL AIRWAYS - SIZE #1 WHITE  
 GUEDEL AIRWAYS - SIZE #2-GREEN  
 GUEDEL AIRWAYS - SIZE #3-YELLOW  
 GUEDEL AIRWAYS - SIZE #4- RED  
 GUEDEL AIRWAYS - SIZE #5-BLUE  
 GUEDEL AIRWAYS - SIZE #6 ORANGE

GUEDEL AIRWAY #1 PORTEX  
 GUEDEL AIRWAY #0 POTEX  
 GUEDEL AIRWAY #00 POTEX  
 GUEDEL AIRWAY #2 PORTEX  
 GUEDEL AIRWAY #3 PORTEX  
 GUEDEL AIRWAY #4 PORTEX

## ENDOTRACHEAL TUBE

### Characteristics

ENDOTRACHEAL TUBE REINFORCED,CUFFED 10.  
 ENDOTRACHEAL TUBE REINFORCED,CUFFED-3.5  
 ENDOTRACHEAL TUBE REINFORCED,CUFFED-4.0  
 ENDOTRACHEAL TUBE REINFORCED,CUFFED-4.5  
 ENDOTRACHEAL TUBE REINFORCED,CUFFED-5.0  
 ENDOTRACHEAL TUBE REINFORCED,CUFFED-5.5  
 ENDOTRACHEAL TUBE REINFORCED,CUFFED-7.0  
 ENDOTRACHEAL TUBE REINFORCED,CUFFED-6.5  
 ENDOTRACHEAL TUBE REINFORCED,CUFFED-8.5  
 ENDOTRACHEAL TUBE REINFORCED,CUFFED-8.0  
 ENDOTRACHEAL TUBE REINFORCED, CUFFED, SIZE  
 ENDOTRACHEAL TUBE REINFORCED,CUFFED-9.5  
 ENDOTRACHEAL TUBE REINFORCED,CUFFED-9.0  
 ENDOTRACHEAL TUBE REINFORCED,CUFFED-6.0  
 ENDOTRACHEAL TUBE STANDARD CUFFED 3.0  
 ENDOTRACHEAL TUBE STANDARD CUFFED 3.5  
 ENDOTRACHEAL TUBE STANDARD CUFFED 4.0  
 ENDOTRACHEAL TUBE STANDARD CUFFED 4.5  
 ENDOTRACHEAL TUBE STANDARD CUFFED 5.0  
 ENDOTRACHEAL TUBE STANDARD CUFFED 5.5



ENDOTRACHEAL TUBE STANDARD,CUFFED 6.0  
 ENDOTRACHEAL TUBE STANDARD,CUFFED 6.5  
 ENDOTRACHEAL TUBE STANDARD,CUFFED 7.0  
 ENDOTRACHEAL TUBE STANDARD,CUFFED 7.5  
 ENDOTRACHEAL TUBE STANDARD,CUFFED 8.0  
 ENDOTRACHEAL TUBE STANDARD,CUFFED 8.5  
 ENDOTRACHEAL TUBE STANDARD,CUFFED 9.5  
 ENDOTRACHEAL TUBE STANDARD,CUFFED 9.0  
 ENDOTRACHEAL TUBE UNCUFFED 3.5  
 ENDOTRACHEAL TUBE UNCUFFED 6.0  
 ENDOTRACHEAL TUBE UNCUFFED 6.5  
 ENDOTRACHEAL TUBE UNCUFFED 7  
 ENDOTRACHEAL TUBE UNCUFFED 8  
 ENDOTRACHEAL TUBE UNCUFFED 8.5  
 ENDOTRACHEAL TUBE UNCUFFED 5.5  
 ENDOTRACHEAL TUBE UNCUFFED 7.5  
 ENDOTRACHEAL TUBE UNCUFFED 4.5 FR  
 ENDOTRACHEAL TUBE UNCUFFED 4 FR  
 ENDOTRACHEAL TUBE UNCUFFED 5FR




E.T. TUBE  
HOLDER



ANESTHETIC  
BREATHING  
SYSTEM



ENTONOX CIRCUITE




HUMIDITY ADAPTOR



BREATHING FILTERS  
& HME FILTERS



R.A.E TUBE  
PERFORMED



## ORAL PERFORMED TUBE

### Characteristics

ORAL PERFORMED TRACHEAL TUBE 8  
ORAL PERFORMED TRACHEAL TUBE 7  
ORAL PERFORMED TRACHEAL TUBE 6  
ORAL PERFORMED TRACHEAL TUBE 5  
ORAL PERFORMED TRACHEAL TUBE 4



## TRACHEOSTOMY TUBE



## LARYNGEAL MASK AIRWAY

### Characteristics

LARYNGEAL MASK AIRWAY - CLASSIC# 2  
LARYNGEAL MASK AIRWAY - CLASSIC# 3  
LARYNGEAL MASK AIRWAY - CLASSIC# 4  
LARYNGEAL MASK AIRWAY - CLASSIC# 5  
LARYNGEAL MASK AIRWAY - FASTR  
LARYNGEAL MASK AIRWAY - FLEXI  
LARYNGEAL MASK (LMA)PROSEAL INTRODUCER

